



myglucohealth[®]

World's Smartest Mobile Diabetes Care Solution



Wireless

**MyGlucoHealth Wireless
transmits glucose
readings via mobile
phone for sharing, review
and analysis.**




MyGlucoHealth
Wireless meter



MyGlucoHealth Transmits Results
Using Bluetooth® Technology



MyGlucoHealth Network Patient Portal
www.myglucohealth.net



Smarter Diabetes Management for a Digital Lifestyle

MyGlucoHealth is an intelligent diabetes management platform providing comprehensive tools for patients, clinicians, families and caregivers to improve lives and lifestyles.

The New MyGlucoHealth® Wireless meter is the world's first integrated Bluetooth-enabled blood glucose meter. With the smallest blood sample size (.3µL), automatic coding and results in just 3 seconds, the MyGlucoHealth Wireless meter is the most technically advanced diabetic care product on the market today.

The MyGlucoHealth Wireless meter along with the MyGlucoHealth Network offers a comprehensive diabetes management platform that gives patients more direct control over their care while providing clear lines of communication with their clinician, family and caregivers.

MyGlucoHealth Wireless works in combination with a real-time online data collection network to upload and report patient blood glucose readings. Using the MyGlucoHealth Patient Portal at www.myglucohealth.net, (or regional sites throughout the world) patients can transmit, review and evaluate their daily readings while communicating those same results automatically to their clinician or caregiver. This allows clinicians to stay up-to-date on patients and intervene immediately if needed.

The MyGlucoHealth Network Portal is an extensive communication and data management web site that provides an easy-to-use interface for automatically collecting test results from the meter. Stored results are wirelessly transmitted by mobile phone or PC to the Portal for posting. This is performed automatically using Bluetooth technology or by connecting a USB cable directly to the meter. The direct upload process eliminates the manual logging of test results, improving accuracy, providing better access to critical data and removing the opportunity for mistakes or "manufactured" entries.

Test, Upload, Manage...Live!



MyGlucoHealth Network Portal

Personal Diabetes Care Just Got Smarter

The MyGlucoHealth Network

- An integrated solution for managing communications between patients and clinicians to effectively manage and treat diabetes.
- Test results are wirelessly transmitted to the Portal using Bluetooth technology.
- The Portal is a private, secure site where patients can review and evaluate results.
- Automated online system for alerts, testing and medication reminders.

Make Better Lifestyle Choices

On the Patient Portal, readings are interpreted and charted in a number of ways to best analyze and evaluate blood glucose data. The Portal is designed to help each patient take "ownership" of their disease, so they can make better lifestyle choices.

Reminders, Messages & Alerts

The MyGlucoHealth Network allows patients to automatically send themselves or their caregivers alerts and reminders about their changing condition or to remind them to test or take medication. Alerts and reminders can be sent by email or SMS text message.

Each patient can create personalized high or low "threshold" settings for blood glucose level. The Network monitors uploaded results and automatically sends an alert to the patient, their family, physician or caregiver when the threshold is exceeded. In addition, patients can send related health information and secure messages directly to their physician and schedule appointments online using the Network Portal.

Patient Reminders

Customizable Targets & Alerts



The MyGlucoHealth Mobile Application is a WAP (wireless application protocol) program that allows users to upload testing seamlessly to the web, as well as review and evaluate their results on their mobile phone.



MyGlucoHealth Wireless

- World's first FDA-cleared Bluetooth® blood glucose meter
- Smallest sample size: 0.3 microliters
- Fastest testing: 3 second results
- Auto-coding
- Alternate site testing
- Plasma-calibrated
- Add blood to test strips for up to one minute
- Easy-to-read display, large clear buttons
- Price of meter and strips comparable to non-connected meters.
- Test Range — 10-600 mg/dL (0.6 - 33.3 mmol/L)
- Power Source — 3 V (Alkaline Battery, 1.5V AAA Size X 2)
- Battery Life — 2,000 Tests

Secure, Accurate, Patient Communication

Data in the MyGlucoHealth Network is secure, accurate, gathered in real-time and accessible to both patients and clinicians 24 hours a day, seven days a week from any Internet-connected computer.

Entra Health Systems received clearance for sale of the MyGlucoHealth Wireless meter in the United States (FDA K081703). The Meter also has CE Mark in Europe (15054) and TGA approval in Australia (ARTG 161202).

Connect to a mobile phone or PC

Transmit readings to the Internet

Share, review and analyze your data

Receive emails, alerts and reminders

Stay connected and in control

Americas

Entra Health Systems LLC
3111 Camino Del Rio North, Suite 101
San Diego, California 92108 U.S.A.
Phone: 1 877 458 2646
1 877 GLUCO 4 ME
1 619 584 6704 (outside the US)
Fax: 413 513 6343
www.myglucohealth.net
www.entrahealthsystems.com
Email: support@myglucohealth.net

Australia and Asia

Entra Health Systems Pty Ltd
12/60 Castlereagh
Sydney, NSW 2000 Australia
Phone: 028 005 4745
Fax: 028 088 7105

Europe, Middle East and Africa

Entra Health Systems Ltd
4A Blake Mews, Kew Gardens
Surrey TW9 3GA United Kingdom
Phone: +44 (0) 113 815 5108